

GOVERNMENT OF TAMILNADU

DIRECTORATE OF TECHNICAL EDUCATION

DIPLOMA IN ENGINEERING I YEAR

SEMESTER SYSTEM
L - SCHEME
2011 - 2012

I SEMESTER

ENGINEERING CHEMISTRY – I PRACTICAL

CURRICULUM DEVELOPMENT CENTER

STATE BOARD OF TECHNICAL EDUCATION & TRAINING, TAMILNADU DIPLOMA IN ENGINEERING - SYLLABUS L-SCHEME

(Implements from the Academic Year 2011-2012 on wards)

Course Name: All Branches of Diploma in Engineering and Technology and Special

Programmes except DMOP, HMCT and Film & TV

Subject Code: 21008

Semester : I Semester

Subject Title : ENGINEERING CHEMISTRY – I PRACTICAL

SCHEME OF INSTRUCTIONS AND EXAMINATION:

No. of Weeks per Semester: 16 Weeks

	Instructions		Examination			
Subject	Hours/ Week	Hours/ Semester	Internal assessment/ Record	Board Examination	Total	Duration
ENGINEERING CHEMISTRY- I PRACTICAL	2	32	25	75	100	3 Hours

OBJECTIVES:

- 1. At the end of the program the student will have knowledge about volumetric analysis in acidimetric, Alkali metric and permanganametric titrations and their applications.
- 2. To give knowledge of estimation of total Hardness, temporary and permanent hardness in the hard water sample.
- 3. To get knowledge about measurement of pH and to calculate Hydrogen ion concentration in a solution.

1. INTERNAL ASSESSMENT/ RECORD: 25 MARKS

2. BOARD EXAMINATION EVALUATION: 75 MARKS

VOLUMETRIC ANALYSIS / EDTA TITRATION

Procedure	5 Marks
Viva-Voce	5 Marks
I-Titration	25 Marks
II Titration	25 Marks
Calculations	3x5 = 15 Marks
Total	75 Marks

DETERMINATION OF PH

Answer for short Question on pH	5 Marks
Viva-voce	5 Marks
Determination of pH (5 Samples)	40 Marks
Calculation of H ⁺ ion concentration	25 Marks
Total	75 Marks

ENGINEERING CHEMISTRY-I PRACTICAL CONTENTS

Intellectual Skills:

- 1. Titrations and Calculation of masses.
- 2. Knowing units for concentration of solutions

Motor Skills:

- 1. Measure the quantities accurately
- 2. Handling the apparatus carefully.

Acidimetery and Alkalimetery:

- Estimation of Sulphuric acid, using a standard solution of oxalic acid and NaOH as Link solution - Phenolphthalein indicator.
 - [Test solution should be made up to 100ml]
- Estimation of sodium hydroxide using a standard solution of sodium carbonate Using sulphuric acid-as link solution -Methyl orange indicator.
 - [Test solution should be made up to 100ml]

- 3. Comparison of strengths of two acid solutions using a link solution of NaOH-Phenolphthalein Indicator.
- Comparison of strengths of two alkaline solutions using an acid (oxalic acid). Phenolphthalein Indicator.

Permanganametry:

- 5. Estimation of Mohr's salt solution using a standard solution of ferrous sulphate and link solution of potassium permanganate
 - [Test solution should be made up to 100ml]
- Estimation of ferrous Sulphate using standard solution of Mohr's salt solution and link solution of potassium permanganate (Test solution should be made up to 100 ml).

Water Analysis:

- 7. Estimation of total hardness of a water sample using EDTA.
- 8. Determination of pH using a pH-meter [for five given samples] and to calculate the hydrogen ion concentration in the solutions.[This question may be given to any two students per batch].

MODEL QUESTION PAPER

MODEL: 1 3 Hours

 Estimate the mass of Sulphuric acid Present in whole of the given solution. You are supplied with a standard solution of oxalic acid of strength 0.098N and an approximately decinormal solution of Sodium hydroxide.

MODEL: 2 3 Hours

2. Calculate the total hardness of the given sample of water. You are given a standard Hard water Solution of 0.01M and an approximately 0.01M EDTA solution.

MODEL: 3 3 Hours

3. Calculate pH of given five samples, using pH meter and calculate the H⁺ ion Concentration of samples. (Any two Students only in a batch).

SCHEME OF VALUATION

Short Procedure (Common to all titration). /

For pH determination question, any two PART A questions -in pH chapter may be asked.

Viva Voce (common to all) - 5 Marks

I. VOLUMETRIC ANALYSIS:-

Titration Value accuracy for I & II----- \pm 0.2 ml -25x2=50

above \pm 0.2 to 0.4 ml - 21 Marks

above \pm 0.4 to 0.6 ml - 17 Marks

above ± 0.6 ml - 5 Marks

Calculations:

Titration I - 5 Marks
Titration II - 5 Marks
Result - 5 Marks

For Arithmetic errors, 25 % marks may be reduced

II. <u>DETERMINATION OF pH VALUES</u>:

Determination of pH for 5 Samples - 5x8=40 Marks

If Accuracy \pm 0.2 - 8 Marks

If Accuracy \pm 0.2 to 0.4 - 6 Marks

If Accuracy \pm 0.4 and above - 4 Marks

Calculation of H^+ ions concentration - 5x5 = 25 Marks

III. EDTA TITRATION

Titration-I - 25 Marks
Titration-II - 25 Marks
Calculation - 15 Marks

For Arithmetic errors, 25 % marks may be reduced

Accuracy on par with ordinary titration (as above)

FIRST YEAR ENGINEERING CHEMISTRY LABORATORY

LIST OF EQUIPMENTS

List of equipments needed for a batch of 30 students in Chemistry Laboratory

Non-Consumable Items:

1.	Indane gas Connection (DBC)	1 no
2.	Exhaust Fan (High capacity)	Sufficient Numbers
3.	Fire Extinguisher	1 no
4.	First Aid Box (Full set)	2 nos
5.	Safety chart	1 no
6.	Chemical Balance	1 no
7.	Fractional weight box	1 no
8.	PH meters	5 nos
9.	Working Table with all accessories	8 nos

Glassware and Other Items:

1.	Burette 50ml	35 nos
2.	Pipette 20ml (with safety bulb)	35 nos
3.	Conical Flask 250ml	35 nos
4.	Funnel 3" (Polythene)	50 nos
5.	Porcelain Tile 6x6"	35 nos
6.	Measuring Cylinder	
	a. 100ml	5 nos
	b. 500 ml	3 nos
7.	Reagent Bottle (White) 250ml	60nos
8.	Reagent Bottle (White) 125ml	100 nos
9.	Reagent Bottle (Amber) 250ml	80 nos

10. Test tubes

a. 15 x 1.5mm	1000 nos
b. 15 x 2.5mm	500 nos
11.Test tube stand	35 nos
12. Test tube holder	35 nos
13. Test tube cleaning brush	35 nos
14. Glass Trough	5 nos
15. Beakers	
a. 1000 ml	5 nos
b. 500 ml	5 nos
c. 250 ml	35 nos
d. 100 ml	5nos
16. Glass Rods 15cm	100 nos
17. Watch Glass 3"	35 nos
18. Wash Bottle (Polythene) 1000ml	35 nos
19. Nickel Spatula	10 nos
20. Kipps Apparatus	1 no
21. Burner Nipple	30 nos
22. Bunsen Burner for gas connection	30 nos
23. Wire Gauge with asbestos center	15 nos
24. Plastic Buckets (15 lts)	10 nos
25. Tripod Stand (Iron)	30 nos
26. Filter Paper Round sheets	1000 nos
27. Burette stand	35 nos
28. Standard flask 100 ml	35 nos
29. Pipette 10ml	5 nos

* * *

List of Apparatus to be supplied for each student in Chemistry lab during FIRST Semester Practical Class/Board Examination in addition to the required Solutions.

1.	Burette 50ml	- 1
2.	Pipette 20ml (with safety bulb)	- 1
3.	Conical Flask 250ml	- 1
4.	Funnel	- 1
5.	Porcelain Tile 6x6"	- 1
6.	Burette stand	- 1
7.	Standard flask 100 ml	- 1
8.	Beakers 250 ml	- 1
9.	Wash Bottle	- 1
		xXx